
INSTRUCTIONS FOR USE

Stride OA - KNEE BRACE

art no 70652

Calf suspension wrap

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EN The wrap prevents rotation of the orthosis during use.

Step 1. Attach self-adhesive Velcro strap on to the orthosis at calf. The round Velcro dot just below the joint, see picture.

Step 2. Attach the wrap with the neoprene side towards skin and pull tight around back of leg.

Step 3. Pull the wrap around the leg and attach.

NOTE! Wrap should be positioned so that top edge of wrap is 1,5 cm above top edge of calf strap 1.

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EN The brace comes with standard 5 degree extension stops. If additional range of motion is needed, use the included stops for 0°, 5°, 10°, 15°, 20°, and 25° extension and 45°, 60°, 75° and 90° flexion.

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EN Mild to severe degree of one compartmental knee arthrosis, gonarthrosis.

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EN Not intended for people with weight >110 kg.

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EN Prescription and treatment recommendations are to be made by a qualified medical practitioner. This also applies to initial fitting. A qualified medical practitioner should also supervise when using this product in combination with other medical products. Follow instruction for best result. Periodically check the condition of the Velcro to ensure they are free from grease and damage. Note that the orthosis is not too tight to avoid stasis. Serious incidents shall be reported to NordiCare and the competent authority.

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EN Salt water, sand, dirt and other elements can cause damage. Thoroughly rinse brace and air dry. Liners and calf wrap should be removed and hand washed using cold water and mild detergent, air dry. Lubricate hinges with dry silicone spray as needed.

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EN Aluminum bar, neoprene knee bag. Lined with suede-like material of polyester. The silicone on the inside for anti-slip.

For more information please consult our product catalogue or visit www.nordicare.se. For further questions, contact your doctor, health practitioner or medical equivalent - or NordiCare: +46 42-35 22 20, info@nordicare.se.

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EN Easy and convenient low profile orthotics used for pain relief and to facilitate and improve function in osteoarthritis of the knee. Unique opportunities to angle the rail and can be positioned for the medial or lateral relief.

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EN Important, before the initial fitting begins:

In the initial fitting, the orthosis should be set in the patient's neutral position (varus or valgus).

1. Place the orthosis on the outside of the leg in a stretched position.
2. Adjust with the Allen key so that the orthosis is fully aligned and follows the patient's individual contour. The desired compartment relief is done in step 5 below.

1. Start:

- 1.1. Loosen straps 1 - 4 and attach each strap back on to itself to facilitate application.

2. Position brace:

- 2.1. Patient to be seated on chair (seated upright with feet planted to floor), bend knee slightly and apply brace to lateral side (outside) of leg.
- 2.2. Align center of hinge with top of kneecap and position slightly behind midline of leg.

3. Fasten straps: in order as indicated on strap ends.

- 3.1. Strap 1 fastens just above calf muscle and is important to maintain brace position.
- 3.2. Strap 2 fastens at mid of calf.
- 3.3. Strap 3 fastens just above knee.
- 3.4. Strap 4 fastens around upper thigh.

All straps should be comfortably tightened to properly secure brace.

Tip - if the patient is experiencing pressure from

the straps: Elastic strap extension may be used for increased comfort:

To relieve medial compartment; open and pull strap 3 out of the buckle. Extend with the elastic extension and thread in to the buckle. To relieve lateral compartment similarly extend strap 4.

4. Check fit:

- 4.1. Once straps are secured, center of hinge should be aligned with middle of kneecap and positioned slightly behind midline of leg.
- 4.2. To ensure proper fit for activity walk a few steps and make sure position is correct.

5. Set the Q angle: Seat patient with knee in 90

degree angle. Use the hex key and adjust the Q-angle an eighth of a turn (1/8):

- 5.1. Clockwise to relieve lateral compartment.
- 5.2. Opposite to relieve medial compartment.

NOTE! Initial fitting and adjustment of the orthosis should be performed by a medically trained person.

6. Further adjustment: Walk and stand to test the orthosis and make corrections as needed. If increased angleation is desired, repeat adjusting the Q-angle by eighth of a turn (1/8) and try again. Excessive correction of the angle is not recommended as in can cause discomfort, increase angulation step by step.

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EN Patient use

Successively increase the time of using the brace, start using the brace for 30 minutes the first day. To obtain best possible effect, wear the brace during activities that strain the knee, such as walking and standing. Arthritis caused pain may take a few weeks to improve.