INFORMATION FOR

ARCUS - CONTRACTURE - HEEL ULCER - AFO

art no 80436-51

Washing instructions

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EN The shell is cleaned with a damp cloth. The padding can be machine washed at 40°C (104°F) with mild detergent. Dried flat while keeping a critical distance from the heat sources.

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EN To relieve/prevent heel ulcer when in bed position. Prevents contracture in the ankle. Counteracts hip rotation.
Contraindications: Unstable ankle fracture. Use caution in case of decreased sensibility in the foot.
Risks: Not intended for walking. Can handle weight during shorter locomotions.

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EN If a stinging sensation, rash or other changes in skin color occur, discontinue use and contact your doctor or other medically trained person. Only medically trained persons are to prescribe and recommend the use of the product. For best results, follow instructions for use. Periodically check the condition of the Velcro to ensure they are free from grease and damage. Use in combination with other medical products should be supervised by qualified medical practitioner. Serious incidents shall be reported to NordiCare and the competent authority.

Material Composition

Plastic rail: 100% Polypropylene.
Padding: 100% Polyurethane (PU) foam with a polyester coating 100% Polyester.
Velcros: 100% Polyester.

For more information, see productleaflet or www.nordicare.se. You can also contact your fitter or NordiCare: +46 42-35 22 20 / DK 7070-1907, info@nordicare.se.

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MD Class I **CE**

 This packaging can be recycled as paper and the bag as plastic.

2022-03-31

NordiCare

Measurements & Size

SE Omkretsmått vaden i centimeter.

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EN The calf circumference measurement is in centimetres.

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



2. Position foot ensuring heel is positioned in recess and toe rest is correct length (see the part re *Adjusting the length/toe rest*).

3. Close fastening on liners ensuring firm but not too tight.

4. Position anti-rotation bar and/or the straps if required (For adjusting the Anti-rotation bar, see the part re *Positioning the Anti-rotation bar*. For adjusting the straps, see the part re *Adjusting dorsiflexion assist straps*.)

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

EN 1. Unfasten the hook/loop flap at heel.

2. Slide the sole off.

3. Pull the toe rest out (There will be some resistance from inner hook fixation). Extend length to accommodate larger foot.

4. Slide sole back into slot. Re-fasten the hook/loop flap at heel. ALWAYS ENSURE SECURE FASTENING

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
EN 1. Unfasten the hook/loop flap at heel **A**. Then open liner **B** as shown.

2. Rotate the bar to desired position. Close liner as shown. ALWAYS ENSURE SECURE FASTENING.

Adjusting dorsiflexion assist straps

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EN 1. Ensure dorsiflexion assist straps are located to front of foot for best leverage. The strap underneath the foot can be easily repositioned with velcro fastening.

2. Tighten straps on both sides to the desired position; dorsiflektion/supination/pronation.